



December 19, 2017

Subject: Cleaning Agents for GE Monitoring Products

Dear Customer,

As preferred cleaning practices change and new cleaning products become available, GE Healthcare continues to test the compatibility of cleaning agents with GE monitoring products. This letter is to provide an overview of the various cleaning agents we have tested and the determination of those tests by GE product.

	B850/B650/B450 Sys	PDM	Central Station V2	Central Station V1	FH Transceiver Sys	ApePro Transmitter Sys	PRN 50M+	Unity Network ID
<b>A = Permitted</b> <b>A<sup>1</sup> = Permitted, can not be used on BX50 touch screens</b> <b>A<sup>2</sup> = Permitted, known to cause battery door to lock after extended use</b> <b>A<sup>3</sup> = Permitted, may cause color shift in soft material after extended use</b> <b>A<sup>4</sup> = Permitted, may cause spotting, staining or corrosion on metal components [example screw, base plate] after extended use</b> <b>A<sup>5</sup> = Permitted, may cause swelling of soft material in battery door after extended use</b> <b>N = Not currently permitted</b>								
Caviwipes Wipes	N	A	A	A	A	A	A	A
Cleanisept Wipes	N	A	A	A	A	A	A <sup>4</sup>	A
Clinell Universal Sanitizing Wipes	N	A	A	A	A	A	A	A
Clinell Sporidical Wipes	A	A <sup>4</sup>	A	A	A	A	A <sup>4</sup>	A
Clorox Hydrogen Peroxide Wipes	A	A <sup>4</sup>	A	A	A <sup>3</sup>	A <sup>3</sup>	A <sup>4</sup>	A <sup>4</sup>
Mikrozyd Sensitive Wipes	N	A	A	A	A	A	A	A
PDI Sani-Cloth AF3 Germicidal Disp Wipes	N	A	A	A	A	A	A <sup>4</sup>	A
PDI Sani-Cloth Bleach Germicidal Disp Wipes	N	A <sup>3</sup>	A	A	A	A <sup>2 3</sup>	A <sup>4</sup>	A <sup>4</sup>
PDI Sani-Cloth HB Germicidal Disp Wipes	N	A	A	A	A	A	A <sup>4</sup>	A <sup>4</sup>
PDI Sani-Cloth Plus Wipes	N	A	A	A	A	A	A	A
PDI Super Sani-Cloth Germicidal Disp Wipes	A	A <sup>3</sup>	A	A	A	A <sup>2</sup>	A	A
PDI Easy Screen Cleaning Wipes	A	A	A	A	A	A	A	A
Oxivir TB Wipes	A	A <sup>4</sup>	A	A	A	A <sup>3</sup>	A	A
Ethanol 70%-96%	A	A	A	N	A	A <sup>5</sup>	A	A
Bleach [dilute to max. of 0.65%][or 0.52%]	A <sup>1</sup>	A <sup>3</sup>	A	A	A	A <sup>2 3</sup>	A <sup>4</sup>	A <sup>4</sup>
Isopropyl alcohol max. 60% by weight	N	A	A	A	A	A	A	A

Note: Identified cleaning agents are permitted for listed products and may not be permitted for accessories. Check accessory user manual for permitted cleaning agents.

DOC2078590

Disinfection Points to Note:

- Always clean before disinfecting.
- Always dilute cleaning and disinfectant agents according to their manufacturer's instructions. Always consider your hospital guidelines as well.
- Use only the permitted substances.
- Your hospital guidelines permitting, all cleaning and disinfection activities can be carried out at the bedside.
- Before drying a surface after wiping it, wait the minimum time required according to the substance manufacturer's instructions.
- Visually ensure that no substance residue remains on the equipment.

Testing of additional GE products, such as the V100 and VC150, and additional cleaning agents are underway and we will update the table above once those tests are complete. We are also in the process of updating the product user manuals to reflect the approved cleaning agents, however for simplicity to our customers we wanted to communicate this information via a consolidated list.

Should you have any questions regarding the contents of this letter, please contact your GE Sales Representative.

Sincerely,



Cory Stahl  
Marketing Manager  
Monitoring Solutions